


Pyroptosis/Inflammation

- Your Expertise
- Our Antibodies
- Accelerated Discovery


www.genetex.com




Adaptor & Sensors

Citation Support Orthogonal Validation Protein Overexpression


Cat. No.	Product Name	Clonality	Reactivity	Applications
GTX116487	AIM2 antibody	Rb pAb	Hu	WB, IHC-P
GTX631645	CARD12 antibody [GT1387]	Ms mAb	Hu	WB
GTX102366	CARD12 antibody [N1], N-term	Rb pAb	Hu	WB
GTX133569	NLRP3 antibody	Rb pAb	Hu, Ms	WB, IHC-P
GTX106313	NLRP3 antibody [C3], C-term	Rb pAb	Hu	WB
GTX102474	TMS1 antibody [N1C3]	Rb pAb	Hu	WB, ICC/IF, IHC-P, IHC-Fr




NLRP3 antibody (GTX133569)



TMS1 antibody [N1C3] (GTX102474)




IL1 beta antibody (GTX74034)




Caspase 4 antibody (GTX134552)

Antibody of Inflammatory Markers


Cat. No.	Product Name	Clonality	Reactivity	Applications
GTX102419	Albumin antibody	Rb pAb	Hu, Ms, Rat, Cat, Dog	WB, ICC/IF, IHC-P, IHC-Fr, IHC
GTX101068	Alpha 1 Microglobulin antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
GTX112815	beta 2 Microglobulin antibody	Rb pAb	Hu, Ms, Rat, Mrt	WB, ICC/IF, IHC-P, IP
GTX10026	C Reactive Protein antibody [C1]	Ms mAb	Hu	WB, ELISA, IHC, Purification, Turbidimetry
GTX101262	C Reactive Protein antibody [N1C3]	Rb pAb	Hu, Rat	WB, IHC-P
GTX64437	Eotaxin antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX52588	Eotaxin antibody [11F7]	Ms mAb	Hu	WB, IHC-P, Neutralizing/Blocking
GTX130513	IL10 antibody	Rb pAb	Hu, Ms	WB, IHC-P, IHC-Fr, ELISA
GTX632359	IL10 antibody [GT5111]	Ms mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr
GTX27842	Insulin antibody	Gp pAb	Hu, Ms	ICC/IF, IHC-P, IHC-Fr
GTX636833	Insulin antibody [HL1393]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX109204	Leptin antibody	Rb pAb	Hu, Ms	WB, IHC-P
GTX26056	Leptin antibody [3G7]	Ms mAb	Hu	WB, ELISA, sELISA
GTX60582	MCP1 / CCL2 antibody [2D8]	Ms mAb	Hu, Ms, Rat, Mk	WB, ICC/IF, IHC-P, FACS, ELISA
GTX81767	MCP1 / CCL2 antibody, C-term	Rb pAb	Hu, Ms	WB, ICC/IF, IHC-P, FACS
GTX108422	MIG antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
GTX135126	Myeloperoxidase antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
GTX75318	Myeloperoxidase antibody [2C7]	Ms mAb	Hu	WB, ICC/IF, IHC-P, IHC-Fr, FACS, ELISA, IHC
GTX47392	PIGF antibody, N-term	Rb pAb	Hu	WB, IHC
GTX102710	TNF Receptor I antibody	Rb pAb	Hu	WB
GTX02834	TNF Receptor II antibody [GT1237]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX54424	TNF Receptor II antibody [80M2]	Ms mAb	Hu, Rat	WB, ICC/IF, IP, Functional Assay, Immunoassay




Albumin antibody (GTX102419)




beta 2 Microglobulin antibody (GTX112815)




C Reactive Protein antibody [N1C3] (GTX101262)




IL10 antibody [GT5111] (GTX632359)




Leptin antibody (GTX109204)



MCP1 / CCL2 antibody [2D8] (GTX60582)




Myeloperoxidase antibody (GTX135126)




MIG antibody (GTX108422)

ELISA kit of Inflammatory Markers

Cat. No.	Product Name	Reactivity
GTX37051	Albumin ELISA kit	Rat
GTX37052	Albumin ELISA kit	MS
GTX37053	Albumin ELISA kit	Hu
GTX37034	beta 2-Microglobulin ELISA Kit	Hu
GTX37033	beta 2-Microglobulin/B2M ELISA Kit	Rat
GTX36938	Prealbumin ELISA Kit	Ms




Albumin ELISA Kit (GTX37053)




beta 2-Microglobulin ELISA Kit (GTX37034)

Proteins of Inflammatory Markers

Cat. No.	Product Name	Reactivity	Applications
GTX00075-pro	Bovine Insulin protein, His and GST tag	Bov	Functional Assay
GTX00242-pro	Human CXCL10 / IP10 protein, His tag (active)	Hu	Functional Assay
GTX00086-pro	Human IL10 protein, His tag (active)	Hu	Functional Assay
GTX00415-pro	Human MCP1 / CCL2 protein, His tag (active)	Hu	Functional Assay
GTX00254-pro	Human PAPP A protein, His tag	Hu	Functional Assay
GTX00234-pro	Human TNF Receptor I protein, His and GST tag (active)	Hu	Functional Assay
GTX00235-pro	Human TNF Receptor II protein, His tag (active)	Hu	Functional Assay
GTX00291-pro	Mouse Insulin protein, His tag (active)	Ms	Functional Assay



Bovine Insulin protein (active) (GTX00075-pro)




Human CXCL10 / IP10 protein (active) (GTX00242-pro)


Effects of IL-10 on IL-1 β production by THP1 cells

Sample (cell supernatant of THP1 cells)	O.D.value	corrected	Concentration of IL-1 β (ng/ml)
Unstimulates	0.151	-0.009	<2.343
Stimulates with LPS (1 μ g/ml)	0.676	0.516	19.9
Stimulates with LPS (1 μ g/ml) and IL-10 (5ng/ml)	0.160	0.000	<2.343
Stimulates with LPS (1 μ g/ml) and IL-10 (10ng/ml)	0.150	-0.010	<2.343


Human IL10 protein (active) (GTX00086-pro)




Human PAPP A protein (active) (GTX00254-pro)




Human MCP1 / CCL2 protein (active)-(GTX00415-pro)




Human TNF Receptor II protein (active) (GTX00235-pro)



Human TNF Receptor I protein (active) (GTX00234-pro)




Mouse Insulin protein (active) (GTX00291-pro)




R Recombinant
 C Citation Support
 K KO/KD Validation
 A Comparable Abs
 M IP/MS Analysis
 O Orthogonal Validation
 P Protein Overexpression

Pyroptosis Signaling


Cat. No.	Product Name	Clonality	Reactivity	Applications
R C A M GTX637906	Caspase 1 antibody [HL2005]	Rec Rb mAb	Hu	WB
C A M GTX111630	Caspase 1 antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
R C A M GTX133447	Caspase 1 (cleaved Asp297) antibody	Rb pAb	Hu	WB
C A M GTX134551	Caspase 1 p10 subunit antibody	Rb pAb	Hu, Ms, Rat	WB, IHC-P
C A M GTX123675	Caspase 1 p10 subunit antibody	Rb pAb	Hu, Ms, Rat	IHC-P
R C A M GTX134552	Caspase 4 antibody	Rb pAb	Hu, Ms, Rat	WB, IHC-P
R C A M GTX113639	Caspase 4 antibody	Rb pAb	Hu, Dog, Pig	WB, IHC-P
R C A M GTX636896	Gasdermin D antibody [HL1430]	Rec Rb mAb	Hu	WB, IHC-P
C A M GTX116840	Gasdermin D antibody [N1N3]	Rb pAb	Hu	WB
C A M R GTX636887	IL1 beta antibody [HL1421]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, ELISA
C A M GTX634188	IL1 beta antibody [GT289]	Ms mAb	Hu	WB, ICC/IF
C A M GTX74034	IL1 beta antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr, ELISA
R GTX637584	IL18 antibody [HL1859]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
R C A M K M A M GTX107678	NFkB p65 antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP, EMSA
K GTX101150	NFkB p100 antibody [C2C3], C-term	Rb pAb	Hu, Ms	WB, IP
R C A M R GTX637436	NFkB p105 antibody [HL1784]	Rec Rb mAb	Hu	WB
K C GTX110585	NFkB p105 antibody	Rb pAb	Hu, Ms, Dog	WB, ICC/IF, IHC-P, ChIP assay, IHC
C GTX104288	P2X7 antibody	Rb pAb	Hu, Ms, Rat	IHC-P, IHC-Fr
GTX124256	P2X7 antibody [C3], C-term	Rb pAb	Ms	WB
GTX130121	Pannexin 1 antibody	Rb pAb	Hu	WB




Gasdermin D antibody [HL1430] (GTX636896) RCAM



NFkB p65 antibody (GTX107678) CKAMROP



IL1 beta antibody [HL1421] (GTX636887) RCAM



Caspase 1 antibody [HL2005] (GTX637906) RCAM



GeneTex

International
USA

Tel: 886.3.6208988
Tel: 1.949.553.1900

Toll-free: 1.877.436.3839

www.genetex.com

PD-P007-230407